





FIBERTECH OPTICA is introducing a range of precision micro optical assemblies featuring optical fibres, microlenses, prisms, GRIN lenses, filters, etc., assembled to high tolerances.

OCT SIDE LOOKING CATHETER

- 500 μm OD
- 8 mm length
- 5 m single mode fiber
- Configurable working distance
- Optional multimode fiber sizes



MINIATURE TRANSMISSION PROBE

- Single pass
- 2 mm path length
- Monolithic assembly
- High throughput
- Less than 3 mm diameter

COLLECTION OPTIC

- High NA optical converter
- Aspheric lens
- Matched with silica 0.22 NA fibers

